



When you need a compact, self-contained cabinet X-ray to provide the sharpest images.

The XPERT[®] 20 Cabinet X-ray System

The most compact,
most powerful
core imaging solution
for science and research
and non-destructive testing.


XPERT 20 System Features:

- AEC based fully automated digital imaging with manual option included.
- Automatic calibration with zero warm up time.
- Ability to save annotated images in multiple formats.
- Seamless interface with pathology and other third party software.
- Contact and 1.3X geometric magnification stages.
- Desktop or mobile platform option for easy transportation.
- Seed counting software.
- Bench-top or tower configurations available.

KUBTEC
S C I E N T I F I C

The XPERT[®] 20

Cabinet X-ray Systems



Specifications

Detector Size	3" x 5" (12 cm x 7 cm)
Spatial Resolution	10.4 lp/mm, contact mode 15.6 lp/mm, mag mode
Energy Range	10-25 kV
Tube Current	up to 1.0 mA
Window Filtration	0.005" beryllium
Power	90-250 VAC, 50/60 Hz, 400 VA
Wireless Connectivity	Standard
Size (W x D x H)	11" x 12" x 20" (28 x 30.5 x 51 cm)
Weight without cart	70 lbs. (32 kg)
Weight with cart	170 lbs. (77 kg)



KUBTEC
S C I E N T I F I C

WWW.KUBTEC.COM    

Specifications subject to change without notice. XPERT, Kubtec and the Kubtec logo are registered trademarks of KUB Technologies, Inc. M1260A-0317